

# Protocol for Silver Staining of Gels



Optimized for Mass Spectrometry and Protein Identification

## GUIDELINES

Silver staining is used for sensitive detection of proteins separated by 1D and 2D SDS PAGE with detection limits from 0.5-5 ng. Many silver staining protocols and commercial staining kits are not compatible with mass spectrometry due to the use of cross-linking reagents. Commercial silver staining kits compatible with mass spectrometry are available from different suppliers, including ProteoSilver Plus (Sigma) and Dodeca Silver Stain (BioRad).

Here we describe a silver staining protocol that has been optimized for mass spectrometric analysis. The protocol results in confident protein identifications and high sequence coverage by MALDI MS and ESI MS due to a high recovery of peptides from the stained gel. The protocol has been tested and documented in many publications (Mortz, E et al, Proteomics 1 (11), 1359-63, 2001).

## MATERIALS

Use high quality chemicals and ultra pure water (18.2 Mohm). Agitate the gels and perform work in a fume hood. Wear gloves at all stages and use clean staining trays with lids to avoid keratin contamination of the gels.

## PROCEDURE

1. Run 1D or 2D gel. Incubate the gel in Fixer (40% ethanol, 10% acetic acid, 50% H<sub>2</sub>O) for 1 hr.
2. Wash the gel in H<sub>2</sub>O for at least 30 min. NB! Overnight washing with several changes of water will remove all acetic acid, reduce background staining and increase sensitivity.
3. Sensitize the gel in 0.02% sodium thiosulfate (0.04 g Na<sub>2</sub>S<sub>2</sub>O<sub>3</sub>, 200 ml H<sub>2</sub>O) for only 1 min. NB! Longer time will decrease peptide recovery from the gel.
4. Wash gel in H<sub>2</sub>O for 3 x 20 sec.
5. Incubate gel for 20 min in 4°C cold 0.1% silver nitrate solution (0.2 g AgNO<sub>3</sub>, 200 ml H<sub>2</sub>O, 0.02% formaldehyde (add 40 µL 35% formaldehyde just before use). NB! Staining is enhanced with cold AgNO<sub>3</sub>.
6. Wash the gel in H<sub>2</sub>O for 3 x 20 sec.
7. Place the gel in a new staining tray. NB! Residual AgNO<sub>3</sub> on the gel surface and staining tray will increase background staining.
8. Wash the gel in H<sub>2</sub>O for 1 min.
9. Develop the gel in 3% sodium carbonate (7.5 g Na<sub>2</sub>CO<sub>3</sub> in 250 ml H<sub>2</sub>O), 0.05% formaldehyde (add 125 µL 35% formaldehyde just before use). NB! Change developer solution immediately when it turns yellow. Terminate when the staining is sufficient.
10. Wash the gel in H<sub>2</sub>O for 20 sec.
11. Terminate staining in 5% acetic acid for 5 min.
12. Leave the gel at 4 °C in 1% acetic for storage. Prior to MS analysis the gel is washed in water for 3 x 10 min to ensure complete removal of acetic acid.

Questions to: [info@alphalyse.com](mailto:info@alphalyse.com)